

Product Data Sheet

Recombinant Human CORO1C Protein

Catalog #	Source	Reactivity	Applications
CST1576	E. coli	H	
Description	Recombinant Human CORO1C Protein is produced by E. coli expression system and the target gene encoding Pro351-Ala470 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CORO1C		
Alternative Names	CRN2; CRNN4; Coronin-1C; Coronin-3; hCRNN4		
Entrez Gene	23603		
SwissProt	Q9ULV4		
Purity			
Chemical Structure	PRKSDLFQDD LYPDTAGPEA ALEAEWFEG KNADPILISL KHGYIPGKNR DLKVVVKNIL DSKPTANKKC DLISIPKKT DTASVQNEAK LDEILKEIKS IKDTICNQDE RISKLEQQMA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com