

Product Data Sheet

Recombinant Human SHCBP1L Protein

Catalog #	Source	Reactivity	Applications
CST1579	E. coli	H	
Description	Recombinant Human SHCBP1L Protein is produced by E. coli expression system and the target gene encoding Thr81-Ala330 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SHCBP1L		
Alternative Names	C1orf14; SHC SH2 domain-binding protein 1-like protein		
Entrez Gene	81626		
SwissProt	Q9BZQ2		
Purity			
Chemical Structure	TGEAAAAAAE EPLLVPPEDE EEAQLPPVC VSRMRGMWRD EKVSLYCDEV LQDCKAEDAD EVMGKYLSEK LKLKDKWLGW WKTNPVSFFV KYEEASIPFV GILVEVTCEP YQDSSSRFKV TVSVAEPFSS NIANIPRDLV DEILEELEHS VPLLEVYPVE GQDTDIHVIA LALEVVRFFY DFLWRDWDDE ESCENYTALE EERINLWCDI QDGTIPGPIA QRFKKTLEKY KN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com