

Recombinant Human SWAP70 Protein

Catalog #	Source	Reactivity	Applications
CST1592	E. coli	H	
Description	Recombinant Human SWAP70 Protein is produced by E. coli expression system and the target gene encoding Val211-Glu580 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SWAP70		
Alternative Names	KIAA0640; Switch-associated protein 70; SWAP-70		
Entrez Gene	23075		
SwissProt	Q9UH65		
Purity			
Chemical Structure	VLKQGYMMKK GHRRKNWTER WFLVKPNIIS YYVSEDLKDK KGDILLDENC CVESLPDKDG KKCLFLVKCF DKTFEISASD KKKKQEWIQA IHSTIHLKL GSPPPHKEAR QRRKELRKKQ LAEQEELERQ MKELQAANES KQQELEAVRK KLEEAASRAA EEEKRLQTQ VELQARFSTE LEREKLIRQQ MEEQVAQKSS ELEQYLQVRV ELEDMYLKLQ EAEDEDERQAR QD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com