

Product Data Sheet

Recombinant Human ACYP2 Protein

Catalog #	Source	Reactivity	Applications
CST1593	E. coli	H	
Description	Recombinant Human ACYP2 Protein is produced by E. coli expression system and the target gene encoding Met1-Tyr99 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ACYP2		
Alternative Names	ACYP; Acylphosphatase-2; Acylphosphatase, muscle type isozyme; Acylphosphate phosphohydrolase 2		
Entrez Gene	98		
SwissProt	P14621		
Purity			
Chemical Structure	MSTAQSLKSV DYEVFGRVQG VCFRMYTEDE ARKIGVVGWV KNTSKGTVTG QVQGPEDKVN SMKSWLSKVG SPSSRIDRTN FSNEKTISKL EYSNFSIRY		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com