

Product Data Sheet

Recombinant Human KLHL41 Protein

Catalog #	Source	Reactivity	Applications
CST1596	E. coli	H	
Description	Recombinant Human KLHL41 Protein is produced by E. coli expression system and the target gene encoding Leu11-Gly380 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	KLHL41		
Alternative Names	KBTBD10; KRP1; Kelch-like protein 41; Kel-like protein 23; Kelch repeat and BTB domain-containing protein 10; Kelch-related protein 1; Sarcosin		
Entrez Gene	10324		
SwissProt	O60662		
Purity			
Chemical Structure	LRLYQSTLLQ DGLKDLLDEK KFDCTLKAG DKSLPCHRLI LSACSPYFRE YFLSEIDEAK KKEVVLDNVD PAILDIIKY LYSASIDLND GNVQDIFALA SRFQIPSVFT VCVSYLQKRL APGNCLAILR LGLLLDCPRL AISAREFVSD RFVQICKEED FMQLSPQELI SVISNDSLNV EKEEAVFEAV MKWVRTDKEN RVKNLSEVFD CIRFRLMTEK YFKDHVEKDD II		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com