

Product Data Sheet

Recombinant Human CINP Protein

Catalog #	Source	Reactivity	Applications
CST1610	E. coli	H	
Description	Recombinant Human CINP Protein is produced by E. coli expression system and the target gene encoding Met1-Leu212 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CINP		
Alternative Names	Cyclin-dependent kinase 2-interacting protein; CDK2-interacting protein		
Entrez Gene	51550		
SwissProt	Q9BW66		
Purity			
Chemical Structure	MEAKTLGTVT PRKPVLVSVA RKIKDNAADW HNLILKWETL NDAGFTTANN IANLKISLLN KDKIELDSSS PASKENEEKV CLEYNEELEK LCEELQATLD GLTKIQVKME KLSSTTKGIC ELENYHYGEE SKRPPLFTW PTTHFYEVSH KLEMYRKEL LLKRTVAKEL AHTGDPDLTL SYLSMWLHQP YVESDSRLHL ESMLLETGHR AL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com