

Product Data Sheet

Recombinant Human FOXL2 Protein

Catalog #	Source	Reactivity	Applications
CST1631	E. coli	H	
Description	Recombinant Human FOXL2 Protein is produced by E. coli expression system and the target gene encoding Met1-Leu376 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FOXL2		
Alternative Names	Forkhead box protein L2		
Entrez Gene	668		
SwissProt	P58012		
Purity			
Chemical Structure	MMASYPEPED AAGALLAPET GRTVKEPEGP PPSPGKGGGG GGGTAPEKPD PAQKPPYSYV ALIAMAIREs AEKRLTSLGI YQYIIAKFPF YEKNKKGWQN SIRHNLSLNE CFIKVPREGG GERKGNWYTL DPACEDMFEEK GNYRRRRRMK RPFRRPPAHF QPGKGLFGAG GAAGGCGVAG AGADGYGYLA PPKYLQSGFL NNSWPLPQPP SPMPYASCQM AAAAAAAAAA AA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com