

Product Data Sheet

Recombinant Human LRRC4C Protein

Catalog #	Source	Reactivity	Applications
CST1635	E. coli	H	
Description	Recombinant Human LRRC4C Protein is produced by E. coli expression system and the target gene encoding Ala361-Asp520 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LRRC4C		
Alternative Names	KIAA1580; NGL1; Leucine-rich repeat-containing protein 4C; Netrin-G1 ligand; NGL-1		
Entrez Gene	57689		
SwissProt	Q9HCJ2		
Purity			
Chemical Structure	ADLNVTEGMA AELKCRASTS LTSVSWITPN GTVMTHGAYK VRIAVLSDGT LNFTNVTVQD TGMVTCMVSN SVGNTTASAT LNVTAATTP FSYFSTVTV TMEPSQDEAR TTDNNVGPTP VVDWETTNVT TSLTPQSTRS TEKTFTIPVT DINSGIPGID		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com