

Product Data Sheet

Recombinant Human DEFB110 Protein

Catalog #	Source	Reactivity	Applications
CST1637	HEK293	H	
Description	Recombinant Human DEFB110 Protein is produced by HEK293 expression system and the target gene encoding Lys20-Gln67 is expressed with C-terminus hFc & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DEFB110		
Alternative Names	DEFB10; DEFB11; DEFB111; Beta-defensin 110; Beta-defensin 10; DEFB-10; Beta-defensin 11; DEFB-11; Beta-defensin 111; Defensin, beta 110; Defensin, beta 111		
Entrez Gene	245913		
SwissProt	Q30KQ9		
Purity			
Chemical Structure	KKKYPEYGSL DLRRECRIGN GQCKNQCHEN EIRIAYCIRP GTHCCLQQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com