

Product Data Sheet

Recombinant Human ADAM8 Protein

Catalog #	Source	Reactivity	Applications
CST1639	HEK293	H	
Description	Recombinant Human ADAM8 Protein is produced by HEK293 expression system and the target gene encoding Cys411-His620 is expressed with C-terminus hFc & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ADAM8		
Alternative Names	MS2; Disintegrin and metalloproteinase domain-containing protein 8; ADAM 8; Cell surface antigen MS2; CD156a		
Entrez Gene	101		
SwissProt	P78325		
Purity			
Chemical Structure	CGNLFVERGE QCDCGPPEDC RNRCCNSTTC QLAEGAQCAH GTCCQECKVK PAGELCRPKK DMCDLEEFCD GRHPECPEDA FQENGTPCSG GYCYNGACPT LAQQCQAFWG PGGQAAEESC FSYDILPGCK ASRYRADMCG VLQCKGGQQP LGRAICIVDV CHALTTEGDT AYEPVPEGTR CGPEKVCWKG RCQDLHVYRS SNCSAQCHNH		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com