

Product Data Sheet

Recombinant Human CDHR4 Protein

Catalog #	Source	Reactivity	Applications
CST1647	E. coli	H	
Description	Recombinant Human CDHR4 Protein is produced by E. coli expression system and the target gene encoding Arg451-Cys550 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CDHR4		
Alternative Names	CDH29; Cadherin-related family member 4; Cadherin-like protein 29		
Entrez Gene	389118		
SwissProt	A6H8M9		
Purity			
Chemical Structure	RVQEDAAPHT LLGSVVGTDMDYPHDNIEYY TSGGPTTFAV DRLSGEVHLL GPLDYEQQRL YRLTVLVIDH GQDQNPNNHLL SGSCITIEV EDVNDHAPEC		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com