

Product Data Sheet

Recombinant Human CLEC1A Protein

Catalog #	Source	Reactivity	Applications
CST1654	HEK293	H	
Description	Recombinant Human CLEC1A Protein is produced by HEK293 expression system and the target gene encoding Thr81-Asp280 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CLEC1A		
Alternative Names	CLEC1; C-type lectin domain family 1 member A; C-type lectin-like receptor 1; CLEC-1		
Entrez Gene	51267		
SwissProt	Q8NC01		
Purity			
Chemical Structure	TGQDTISQME ERLGNTSQEL QSLQVQNIKL AGSLQHVAEK LCRELYNKAG AHRCSPCTEQ WKWHGDNQCYQ FYKDSKSWED CKYFCLSENS TMLKINKQED LEFAASQSYS EFFYSYWTGL LRPDSGKAWL WMDGTPFTSE LFHIIIDVTS PRSRDCVAIL NGMIFSKDCK ELKRCVCERR AGMVKPESLH VPPETLGEGD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com