

Recombinant Human LRFN1 Protein

Catalog #	Source	Reactivity	Applications
CST1669	E. coli	H	
Description	Recombinant Human LRFN1 Protein is produced by E. coli expression system and the target gene encoding Arg221-Asn420 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LRFN1		
Alternative Names	KIAA1484; SALM2; Leucine-rich repeat and fibronectin type III domain-containing protein 1; Synaptic adhesion-like molecule 2		
Entrez Gene	57622		
SwissProt	Q9P244		
Purity			
Chemical Structure	RLHKLPPDGL FLRSQGTGPK PPTPLTVSFG GNPLHCNCEL LWLRRLTRED DLETCATPEH LTDRYFWSIP EEEFLCEPPL ITRQAGGRAL VVEGQAVSLR CRAVGDPPEPV VHWVAPDGRL LGNSSRTRVR GDGTLDVTIT TLRDSGTFTC IASNAAGEAT APVEVCVVPL PLMAPPPAAP PPLTEPGSSD IATPGRPGAN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com