

Product Data Sheet

Recombinant Human IGFBPL1 Protein

Catalog #	Source	Reactivity	Applications
CST1676	HEK293	H	
Description	Recombinant Human IGFBPL1 Protein is produced by HEK293 expression system and the target gene encoding Cys41-His270 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	IGFBPL1		
Alternative Names	IGFBPRP4; Insulin-like growth factor-binding protein-like 1; IGFBP-related protein 10; Insulin-like growth factor-binding-related protein 4; IGFBP-rP4		
Entrez Gene	347252		
SwissProt	Q8WX77		
Purity			
Chemical Structure	CRPEGCPAPA PCPAPGISAL DECGCCARCL GAEGASCGGR AGGRCGPGLV CASQAAGAAP EGTGLCVCAQ RGTVCSDGR SYPSVCALRL RARHTPRAHP GHLHKARDGP CEFAPVVVVP PRSVHNVGTGA QVGLSCEVRA VPTPVITWRK VTKSPEGTQA LEELPGDHVN IAVQVRGGPS DHEATAWILI NPLRKEDEGV YQCHAANMVG EAESHSTVTV LDLSKYRSFH		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com