

Product Data Sheet

Recombinant Human IL17D Protein

Catalog #	Source	Reactivity	Applications
CST1677	HEK293	H	
Description	Recombinant Human IL17D Protein is produced by HEK293 expression system and the target gene encoding Ala16-Pro202 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	IL17D		
Alternative Names	IL27; Interleukin-17D; IL-17D; Interleukin-27; IL-27		
Entrez Gene	53342		
SwissProt	Q8TAD2		
Purity			
Chemical Structure	AGAPRAGRPP ARPRGCADRP EELLEQLYGR LAAGVLSAFH HTLQLGPREQ ARNASCPAGG RPADRRFRPP TNLRVSPWA YRISYDPARY PRYLPEAYCL CRGCLTGLFG EEDVRFRSAP VYMPTVVLRR TPACAGGRSV YTEAYVTIPV GCTCVPEPEK DADSINSSID KQGAKLLGP NDAPAGP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com