

Product Data Sheet

Recombinant Human RTN4RL2 Protein

Catalog #	Source	Reactivity	Applications
CST1686	E. coli	H	
Description	Recombinant Human RTN4RL2 Protein is produced by E. coli expression system and the target gene encoding Pro271-Arg380 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RTN4RL2		
Alternative Names	NGRH1; NGRL3; Reticulon-4 receptor-like 2; Nogo receptor-like 3; Nogo-66 receptor homolog 1; Nogo-66 receptor-related protein 2; NgR2		
Entrez Gene	349667		
SwissProt	Q86UN3		
Purity			
Chemical Structure	PLWAWFQRRAR VSSSDVTCAT PPERQGRDLR ALREADFQAC PPAAPTRPGS RARGNSSSNH LYGVAEAGAP PADPSTLYRD LPAEDSRGRQ GGDAPTEDDY WGGYGGEDQR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com