

Product Data Sheet

Recombinant Human ADAMTS4 Protein

Catalog #	Source	Reactivity	Applications
CST1688	E. coli	H	
Description	Recombinant Human ADAMTS4 Protein is produced by E. coli expression system and the target gene encoding Lys441-Asn570 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ADAMTS4		
Alternative Names	KIAA0688; A disintegrin and metalloproteinase with thrombospondin motifs 4; ADAM-TS 4; ADAM-TS4; ADAMTS-4; ADMP-1; Aggrecanase-1		
Entrez Gene	9507		
SwissProt	O75173		
Purity			
Chemical Structure	KDYDADRQCQ LTFGPDSRHC PQLPPPCAAL WCSGHLNGHA MCQTKHSPWA DGTPCGPAQA CMGGRCLHMD QLQDFNIPQA GGWGPWGPWG DCSRTC GGGV QFSSRDCTR VP RN G G KYCE GR TR FR SCN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com