

Product Data Sheet

Recombinant Human TAC3 Protein

Catalog #	Source	Reactivity	Applications
CST1694	HEK293	H	
Description	Recombinant Human TAC3 Protein is produced by HEK293 expression system and the target gene encoding Ala21-Ala120 is expressed with N-terminus 9×His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TAC3		
Alternative Names	NKNB; Tachykinin-3; ZNEUROK1		
Entrez Gene	6866		
SwissProt	Q9UHF0		
Purity			
Chemical Structure	AVCKEPQEEV VPGGGRSKRD PDLYQLLQRL FKSHSSLEGL LKALSQASTD PKESTSPEKR DMHDFVGLM GKRSVQPDSP TDVNQENVPS FGILKYPPRA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com