

Product Data Sheet

Recombinant Human SRPK2 Protein

Catalog #	Source	Reactivity	Applications
CST1709	E. coli	H	
Description	Recombinant Human SRPK2 Protein is produced by E. coli expression system and the target gene encoding Pro321-Gly460 is expressed with N-terminus MBP & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SRPK2		
Alternative Names	SRSF protein kinase 2; SFRS protein kinase 2; Serine/arginine-rich protein-specific kinase 2; SR-protein-specific kinase 2		
Entrez Gene	6733		
SwissProt	P78362		
Purity			
Chemical Structure	PEVKLKTTGL EEAAEAETAK DNGEAEDQEE KEDAEKENIE KEDDDVDQEL ANIDPTWIES PKTNGHIENG PFSLEQQQLDD EDDDEEDCPN PEEYNLDEPN AESDITYSSS YEQFNGELPN GRHKIPESQF PEFSTSLFSG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com