

Product Data Sheet

Recombinant Human GCET1 Protein

Catalog #	Source	Reactivity	Applications
CST1710	E. coli	H	
Description	Recombinant Human GCET1 Protein is produced by E. coli expression system and the target gene encoding Phe331-Lys410 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SERPINA9		
Alternative Names	GCET1; SERPINA11; Serpin A9; Centerin; Germinal center B-cell-expressed transcript 1 protein		
Entrez Gene	327657		
SwissProt	Q86WD7		
Purity			
Chemical Structure	FDKNADFSGI AKRDSLQVSK ATHKAVLDVS EEGTEATAAT TTKFIVRSKD GPSYFTVSFN RTFLMMITNK ATDGILFLGK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com