

Product Data Sheet

Recombinant Human TCL1B Protein

Catalog #	Source	Reactivity	Applications
CST1712	E. coli	H	
Description	Recombinant Human TCL1B Protein is produced by E. coli expression system and the target gene encoding Val11-Leu120 is expressed with N-terminus MBP & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TCL1B		
Alternative Names	TCL1; T-cell leukemia/lymphoma protein 1B; Oncogene TCL-1B; Oncogene TCL1B; SYN-1; Syncytiotrophoblast-specific protein; TCL1/MTCP1-like protein 1		
Entrez Gene	9623		
SwissProt	O95988		
Purity			
Chemical Structure	VPPGRLWIQR PGIYEDEEGR TWVTVVVRFN PSRREWARAS QGSRYEPSIT VHLWQMAVHT RELLSSGQMP FSQLPVWQL YPGRKYRAAD SSFWEIADHG QIDSMEQLVL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com