

Product Data Sheet

Recombinant Human SBSPON Protein

Catalog #	Source	Reactivity	Applications
CST1713	E. coli	H	
Description	Recombinant Human SBSPON Protein is produced by E. coli expression system and the target gene encoding Thr151-Val250 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SBSPON		
Alternative Names	C8orf84; RPESP; Somatomedin-B and thrombospondin type-1 domain-containing protein; RPE-spondin		
Entrez Gene	157869		
SwissProt	Q8IVN8		
Purity			
Chemical Structure	TRQATSPHWS THTEDAGYCM EFKTESLTPH CALENWPLTR WMQYLREGYT VCVDCQPPAM NSVSLRCSGD GLDSDGNQTL HWQAIGNPRC QGTWKKVRRV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com