

Product Data Sheet

Recombinant Human LILRA5 Protein

Catalog #	Source	Reactivity	Applications
CST1721	E. coli	H	
Description	Recombinant Human LILRA5 Protein is produced by E. coli expression system and the target gene encoding Ala51-Asp260 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LILRA5		
Alternative Names	ILT11; LILRB7; LIR9; Leukocyte immunoglobulin-like receptor subfamily A member 5; CD85 antigen-like family member F; Immunoglobulin-like transcript 11; ILT-11; Leukocyte immunoglobulin-like receptor 9; LIR-9; CD85f		
Entrez Gene	353514		
SwissProt	A6NI73		
Purity			
Chemical Structure	AEPGSVISRG NSVTIRCQGT LEAQEYRLVK EGSPEPWDTQ NPLEPKNKAR FSIPSMTEHH AGRYRCYYYS PAGWSEPSDP LELVVTGFYN KPTLSALPSP VVTSGENVTL QCGSRLRFDR FILTEEGDHK LSWTLDSQLT PSGQFQALFP VGPVTPSHRW MLRCYGSRRH ILQVWSEPSD LLEIPVSGAA DNLSPSQNKS DSGTASHLQD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com