

Product Data Sheet

Recombinant Human PAR6 Protein

Catalog #	Source	Reactivity	Applications
CST1730	E. coli	H	
Description	Recombinant Human PAR6 Protein is produced by E. coli expression system and the target gene encoding Gly261-Asp340 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	PAR6A		
Alternative Names	PAR6A; Partitioning defective 6 homolog alpha; PAR-6; PAR-6 alpha; PAR-6A; PAR6C; Tax interaction protein 40; TIP-40		
Entrez Gene	50855		
SwissProt	Q9NPB6		
Purity			
Chemical Structure	GRLTGPPSAG PGPAEPDSDD DSSDLVIENR QPPSSNGLSQ GPPCWDLHPG CRHPGTRSSL PSLDDQEQAS SGWGSRIRGD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com