

Product Data Sheet

Recombinant Human MFRP Protein

Catalog #	Source	Reactivity	Applications
CST1734	E. coli	H	
Description	Recombinant Human MFRP Protein is produced by E. coli expression system and the target gene encoding Gln131-Asp290 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MFRP		
Alternative Names	Membrane frizzled-related protein; Membrane-type frizzled-related protein		
Entrez Gene	83552		
SwissProt	Q9BY79		
Purity			
Chemical Structure	QQESGVSPSP QSTCGLLSG PRGFFSSPNY PDPYPPNTHC VWHIQVATDH AIQLKIEALS IESVASCLFD RLELSPEPEG PLLRVCGRVP PPTLNTNASH LLVVFVSDSS VEGFGFHAWY QAMAPGRGSC AHDEFRCDDL ICLLPDSVCD GFANCADGSD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com