

Product Data Sheet

Recombinant Human HAVCR2 Protein

Catalog #	Source	Reactivity	Applications
CST1739	E. coli	H	
Description	Recombinant Human HAVCR2 Protein is produced by E. coli expression system and the target gene encoding Gly31-Arg200 is expressed with N-terminus GST tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	HAVCR2		
Alternative Names	TIM3; TIMD3; Hepatitis A virus cellular receptor 2; HAVcr-2; T-cell immunoglobulin and mucin domain-containing protein 3; TIMD-3; T-cell immunoglobulin mucin receptor 3; TIM-3; T-cell membrane protein 3		
Entrez Gene	84868		
SwissProt	Q8TDQ0		
Purity			
Chemical Structure	GQNAYLPCFY TPAAPGNLVP VCWKGKACPV FECGNVVLRT DERDVNYWTS RYWLNQDFRK GDVSLTIENV TLADSGIYCC RIQIPGIMND EKFNKLKVIK PAKVTPAPTR QRDFTAAFPR MLTTRGHGPA ETQTLGSLPD INLTQISTLA NELRDSRLAN DLRDSGATIR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com