

Product Data Sheet

Recombinant Human BIK Protein

Catalog #	Source	Reactivity	Applications
CST1747	E. coli	H	
Description	Recombinant Human BIK Protein is produced by E. coli expression system and the target gene encoding Ile11-Gly110 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	BIK		
Alternative Names	NBK; Bcl-2-interacting killer; Apoptosis inducer NBK; BIP1; BP4		
Entrez Gene	638		
SwissProt	Q13323		
Purity			
Chemical Structure	ILMETLLYEQ LLEPPTMEVL GMTDSEEDLD PMEDFDSLEC MEGSDALALR LACIGDEMDV SLRAPRLAQL SEVAMHSLGL AFIYDQTEDI RDVLRSFMDG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com