

Product Data Sheet

Recombinant Human ADGRF4 Protein

Catalog #	Source	Reactivity	Applications
CST1752	HEK293	H	
Description	Recombinant Human ADGRF4 Protein is produced by HEK293 expression system and the target gene encoding Lys31-Leu280 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GPR115		
Alternative Names	PGR18; Probable G-protein coupled receptor 115; G-protein coupled receptor PGR18		
Entrez Gene	221393		
SwissProt	Q8IZF3		
Purity			
Chemical Structure	KAGDKLQSPE GKPKTGRIQE KCEGPCISSS NCSQPCAKDF HGEIGFTCNQ KKWQKSAETC TSLSVEKLFK DSTGASRLSV AAPSIPLHIL DFRAPETIES VAQGIRKNCP FDYACITDMV KSSETTSNGI AFIVELLKNI STDLSDNVTR EKMKSYS EVA NHILD TAAIS NWAFIPNKNA SSDLLQSVNL FARQLHIHNN SENIVNELFI QTKGFHINHN TSEKSLNFSM SM		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com