

Product Data Sheet

Recombinant Human CREG1 Protein

Catalog #	Source	Reactivity	Applications
CST1754	HEK293	H	
Description	Recombinant Human CREG1 Protein is produced by HEK293 expression system and the target gene encoding Arg32-Gln220 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CREG1		
Alternative Names	CREG; Protein CREG1; Cellular repressor of E1A-stimulated genes 1		
Entrez Gene	8804		
SwissProt	O75629		
Purity			
Chemical Structure	RGGRDHGDWD EASRLPPLPP REDAARVARF VTHVSDWGAL ATISTLEAVR GRPFADVLSL SDGPPGAGSG VPHYLYSPLQ LSVSNLQENP YATLTMTLAQ TNFCKKHGFD PQSPLCVHIM LSGTVTKVNE TEMDIAKHSL FIRHPEMKTW PSSHNWFFAK LNITNIWVLD YFGGPKIVTP EEYYNVTVQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com