

Product Data Sheet

Recombinant Human CRISP2 Protein

Catalog #	Source	Reactivity	Applications
CST1757	E. coli	H	
Description	Recombinant Human CRISP2 Protein is produced by E. coli expression system and the target gene encoding Lys22-Tyr243 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CRISP2		
Alternative Names	GAPDL5; TPX1; Cysteine-rich secretory protein 2; CRISP-2; Cancer/testis antigen 36; CT36; Testis-specific protein TPX-1		
Entrez Gene	7180		
SwissProt	P16562		
Purity			
Chemical Structure	KDPAFTALLT TQLQVQREIV NKHNELRKAV SPPASNMLKM EWSREVTNA QRWANKCTLQ HSDPEDRKTS TRCGENLYMS SDPTSWSSAI QSWYDEILDF VYGVGPKSPN AVVGHYTQLV WYSTYQVGCG IAYCPNQDSL KYYYVCQYCP AGNNMNRKNT PYQQGTPCAG CPDDCDKGLC TNSCQYQDLL SNCDLKNNTA GCEHELLKEK CKATCLCENK IY		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com