

Product Data Sheet

Recombinant Human RANGRF Protein

Catalog #	Source	Reactivity	Applications
CST1762	E. coli	H	
Description	Recombinant Human RANGRF Protein is produced by E. coli expression system and the target gene encoding Val31-Pro160 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RANGRF		
Alternative Names	MOG1; RANGNRF; Ran guanine nucleotide release factor; RanGNRF; Ran-binding protein MOG1		
Entrez Gene	29098		
SwissProt	Q9HD47		
Purity			
Chemical Structure	VPDNQEVFCH PVTDAQSLIVE LLELQAHVRG EAAARYHFED VGGVQGARAV HVESVQPLSL ENLALRGRCQ EAWVLSGKQQ IAKENQQVAK DVTLHQALLR LPQYQTDLLL TFNQPPPDNR SSLGPENLSP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com