

Product Data Sheet

Recombinant Human CHRDL2 Protein

Catalog #	Source	Reactivity	Applications
CST1765	HEK293	H	
Description	Recombinant Human CHRDL2 Protein is produced by HEK293 expression system and the target gene encoding Arg26-Thr429 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CHRDL2		
Alternative Names	BNF1; CHL2; Chordin-like protein 2; Breast tumor novel factor 1; BNF-1; Chordin-related protein 2		
Entrez Gene	25884		
SwissProt	Q6WN34		
Purity			
Chemical Structure	RARPDMFCLF HGKRYSPGES WHPYLEPQGL MYCLRCTCSE GAHVSCYRLH CPPVHCPQPV TEPQQCCPKC VEPHTPSGLR APPKSCQHNG TMYQHGEIFS AHELFP SRLP NQCVLC SCTE GQIYCGLTTC PEPGCPAPLP LPDSCCQACK DEASEQSDEE DSVQSLHGVR HPQDPCSSDA GRKRGP GTPA PTGLSAPLSF IPRHFRPKGA GSTTVKIVLK EKHKKACVHG GK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com