

Product Data Sheet

Recombinant Human DEFB127 Protein

Catalog #	Source	Reactivity	Applications
CST1766	E. coli	H	
Description	Recombinant Human DEFB127 Protein is produced by E. coli expression system and the target gene encoding Gly31-Asn90 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DEFB127		
Alternative Names	C20orf73; DEFB27; Beta-defensin 127; Beta-defensin 27; DEFB-27; Defensin, beta 127		
Entrez Gene	140850		
SwissProt	Q9H1M4		
Purity			
Chemical Structure	GHCRKICRVN EVPEALCENG RYCCLNIKEL EACKKITKPP RPKPATLALT LQDYVTIEN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com