

Product Data Sheet

Recombinant Human LRRTM1 Protein

Catalog #	Source	Reactivity	Applications
CST1773	E. coli	H	
Description	Recombinant Human LRRTM1 Protein is produced by E. coli expression system and the target gene encoding Met271-Glu420 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LRRTM1		
Alternative Names	Leucine-rich repeat transmembrane neuronal protein 1		
Entrez Gene	347730		
SwissProt	Q86UE6		
Purity			
Chemical Structure	MEPHVFETVP HLQSLQLDSN RLTYIEPRIL NSWKSLTSIT LAGNLWDCGR NVCALASWLN NFQGRYDGNL QCASPEYAQG EDVLDAVYAF HLCEDGAAPT SGHLLSAVTN RSDLGPPASS ATTLADGGEG QHDGTFEPAT VALPGGEHAE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com