

Product Data Sheet

Recombinant Human SLURP2 Protein

Catalog #	Source	Reactivity	Applications
CST1780	E. coli	H	
Description	Recombinant Human SLURP2 Protein is produced by E. coli expression system and the target gene encoding Ile23-Asp97 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SLURP2		
Alternative Names	Secreted Ly-6/uPAR domain-containing protein 2; Secreted LY6/PLAUR domain-containing protein 2; Secreted Ly-6/uPAR-related protein 2		
Entrez Gene	432355		
SwissProt	P0DP57		
Purity			
Chemical Structure	IWCHQCTGFG GCSHGSRCLR DSTHCVTTAT RVLSNTEDLP LVTKMCHIGC PDIPSLGLGP YVSIACCQTS LCNHD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com