

Product Data Sheet

Recombinant Human TFA1 Protein

Catalog #	Source	Reactivity	Applications
CST1781	E. coli	H	
Description	Recombinant Human TFA1 Protein is produced by E. coli expression system and the target gene encoding Phe31-His130 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FAM19A1		
Alternative Names	TFA1; Protein FAM19A1; Chemokine-like protein TFA-1		
Entrez Gene	407738		
SwissProt	Q7Z5A9		
Purity			
Chemical Structure	FQQHHLHRPE GGTCEVIAAH RCCNKNRIEE RSQTVKCSCL PGKVAGTTRN RPSCVDASIV IGKWWCEMEP CLEGEECKTL PDNSGWMCAT GNKIKTTRIH		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com