

Product Data Sheet

Recombinant Human MDGA1 Protein

Catalog #	Source	Reactivity	Applications
CST1790	E. coli	H	
Description	Recombinant Human MDGA1 Protein is produced by E. coli expression system and the target gene encoding Ile341-Gln530 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MDGA1		
Alternative Names	MAMDC3; MAM domain-containing glycosylphosphatidylinositol anchor protein 1; GPI and MAM protein; GPIM; Glycosylphosphatidylinositol-MAM; MAM domain-containing protein 3		
Entrez Gene	266727		
SwissProt	Q8NFP4		
Purity			
Chemical Structure	IKESENIQLG QDLKLSCHVD AVPQEKVTYQ WFKNGKPARM SKRLLVTRND PELPAVTSSL ELIDLHFSYD GTYLCMASFP GAPVPDLSVE VNISSETVPP TISVPKGRAV VTVREGSPA LQCEVRGKPR PPVLWSRVDK EAALLPSGLP LEETPDGKLR LERVSRDMSG TYRCQTARYN GFNVRPREAQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com