

Product Data Sheet

Recombinant Human NBL1 Protein

Catalog #	Source	Reactivity	Applications
CST1792	HEK293	H	
Description	Recombinant Human NBL1 Protein is produced by HEK293 expression system and the target gene encoding Ala17-Asp181 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	NBL1		
Alternative Names	DAN; DAND1; Neuroblastoma suppressor of tumorigenicity 1; DAN domain family member 1; Protein N03; Zinc finger protein DAN		
Entrez Gene	4681		
SwissProt	P41271		
Purity			
Chemical Structure	AAPPPINKLA LFPDKSAWCE AKNITQIVGH SGCEAKSIQN RACLGQCFSY SVPNTFPQST ESLVHCDSCM PAQSMWEIVT LECPGHEEVP RVDKLVEKIL HCSCQACGKE PSHEGLSVYV QGEDGPGSQP GTHPHPHPHP HPGGQTPEPE DPPGAPHTEE EGAED		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com