

Product Data Sheet

Recombinant Human UNC5C Protein

Catalog #	Source	Reactivity	Applications
CST1793	E. coli	H	
Description	Recombinant Human UNC5C Protein is produced by E. coli expression system and the target gene encoding Glu51-Tyr380 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	UNC5C		
Alternative Names	UNC5H3; Netrin receptor UNC5C; Protein unc-5 homolog 3; Protein unc-5 homolog C		
Entrez Gene	8633		
SwissProt	O95185		
Purity			
Chemical Structure	ETFPSPPEP LPHFLIEPEE AYIVKNKPVN LYCKASPATQ IYFKCNSEWV HQKDHIVDER VDETSGLIVR EVSIEISRQQ VEELFGPEDY WCQCVAWSSA GTTKSRKAYV RIAYLRKTFE QEPLGKEVSL EQEVLLQCRP PEGIPVAEVE WLKNEDIIDP VEDRNFYITI DHNLIQKAR LSDTANYTCV AKNIVAKRKS TTATVIVYVN GGWSTWTEWS VCNSRCGRGY QK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com