

Product Data Sheet

Recombinant Human GAS6 Protein

Catalog #	Source	Reactivity	Applications
CST1796	HEK293	H	
Description	Recombinant Human GAS6 Protein is produced by HEK293 expression system and the target gene encoding Phe41-Cys670 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GAS6		
Alternative Names	AXLLG; Growth arrest-specific protein 6; GAS-6; AXL receptor tyrosine kinase ligand		
Entrez Gene	2621		
SwissProt	Q14393		
Purity			
Chemical Structure	FLRPRQRRAF QVFEEAKQGH LERECVEELC SREEAREVFE NDPETDYFYP RYLDCKINKYG SPYTKNSGFA TCVQNLDPQC TPNPCDRKGT QACQDLMGNF FCLCKAGWGG RLCDKDVNEC SQENGGCLQI CHNKPFSFHC SCHSGFELSS DGRTCQDIDE CADSEACGEA RCKNLPGSYS CLCDEGFAYS SQEKACRDVD ECLQGRCEQV CVNSPGSYTC HCDGRGGLKL SQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com