

Product Data Sheet

Recombinant Human TNFSF14 Protein

Catalog #	Source	Reactivity	Applications
CST1799	HEK293	H	
Description	Recombinant Human TNFSF14 Protein is produced by HEK293 expression system and the target gene encoding Glu81-Val240 is expressed with N-terminus hFc & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TNFSF14		
Alternative Names	HVEM; LIGHT; Tumor necrosis factor ligand superfamily member 14; Herpes virus entry mediator ligand; HVEM-L; Herpesvirus entry mediator ligand; CD258		
Entrez Gene	8740		
SwissProt	O43557		
Purity			
Chemical Structure	EQLIQERRSH EVNPA AHLTG ANSSLTGSGG PLLWETQLGL AFLRGLSYHD GALVVTKAGY YYISKVQLG GVGCPGLGLAS TITHGLYKRT PRYP EELELL VSQQSPCGRA TSSSRVWWDS SFLGGVVHLE AGEKV VVRVL DERLVRLRDG TRSYFGAFMV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com