

Product Data Sheet

Recombinant Human HECTD4 Protein

Catalog #	Source	Reactivity	Applications
CST1802	E. coli	H	
Description	Recombinant Human HECTD4 Protein is produced by E. coli expression system and the target gene encoding Ser3071-Glu3210 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	HECTD4		
Alternative Names	C12orf51; KIAA0614; Probable E3 ubiquitin-protein ligase HECTD4; HECT domain-containing protein 4		
Entrez Gene	283450		
SwissProt	Q9Y4D8		
Purity			
Chemical Structure	SICSSQGISQ TVSDLSVDPL PAGLELPIPP GLLEPHAVSS QESLDISLCS TGSLSGLGSL GEPLDNAETA SVSDMGSMYT VTSLDNQPLA ARPIKGFVAV EIRSRKIEK IRASLFNNND LIGLSSLDGE DELMEMSTEE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com