

Product Data Sheet

Recombinant Human RNF150 Protein

Catalog #	Source	Reactivity	Applications
CST1803	E. coli	H	
Description	Recombinant Human RNF150 Protein is produced by E. coli expression system and the target gene encoding His301-Asp430 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RNF150		
Alternative Names	KIAA1214; RING finger protein 150		
Entrez Gene	57484		
SwissProt	Q9ULK6		
Purity			
Chemical Structure	HKSCVDPWLL DHRTCPMCKM NILKALGIPP NADCMDDLPT DFEGSLGGPP TNQITGASDT TVNESSVTLD PAVRTVGALQ VVQDTPPIQ EGDVIFTTNS EQEPAVSSDS DISLIMAMEV GLSDVELSTD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com