

Product Data Sheet

Recombinant Human RPRD1B Protein

Catalog #	Source	Reactivity	Applications
CST1805	E. coli	H	
Description	Recombinant Human RPRD1B Protein is produced by E. coli expression system and the target gene encoding Glu141-Thr310 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RPRD1B		
Alternative Names	C20orf77; CREPT; Regulation of nuclear pre-mRNA domain-containing protein 1B; Cell cycle-related and expression-elevated protein in tumor		
Entrez Gene	58490		
SwissProt	Q9NQG5		
Purity			
Chemical Structure	EEKSLKRTF QQIQEEEDDD YPGSYSPQDP SAGPLLTEEL IKALQDLENA ASGDATVRQK IASLPQEVQD VSLLEKITDK EAAERLSKTV DEACLLLA EY NGRLLAAELED RRQLARMLVE YTQNQKDVLS EKEKKLEEYK QKLARVTQVR KELKSHIQSL PDL SLLPNVT		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com