

Product Data Sheet

Recombinant Human TRIM31 Protein

Catalog #	Source	Reactivity	Applications
CST1808	E. coli	H	
Description	Recombinant Human TRIM31 Protein is produced by E. coli expression system and the target gene encoding Lys121-Ser390 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TRIM31		
Alternative Names	C6orf13; RNF; E3 ubiquitin-protein ligase TRIM31; Tripartite motif-containing protein 31		
Entrez Gene	11074		
SwissProt	Q9BZY9		
Purity			
Chemical Structure	KDHKSHNVSL IEEAAQNYQG QIQEQIQVLQ QKEKETVQVK AQGVHRVDVF TDQVEHEKQR ILTEFELLHQ VLEEEKNFLL SRIYWLGHG TEAGKHVVAS TEPQLNDLKK LVDSLTKQN MPPRQLLEDI KVVLCRSEEF QFLNPTPVPL ELEKKLSEAK SRHDSITGSL KFKDQLQAD RKKDENRFFK SMNKNDMKSW GLLQKNNHKM NKTSEPGSSS AGGRTTSGPP NH		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com