

Product Data Sheet

Recombinant Human UPF3B Protein

Catalog #	Source	Reactivity	Applications
CST1809	E. coli	H	
Description	Recombinant Human UPF3B Protein is produced by E. coli expression system and the target gene encoding Asp41-Glu420 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	UPF3B		
Alternative Names	RENT3B; UPF3X; Regulator of nonsense transcripts 3B; Nonsense mRNA reducing factor 3B; Up-frameshift suppressor 3 homolog B; hUpf3B; Up-frameshift suppressor 3 homolog on chromosome X; hUpf3p-X		
Entrez Gene	65109		
SwissProt	Q9BZI7		
Purity			
Chemical Structure	DRNKEKKEAL SKVVIRRLPP TLTKEQLQEH LQPMPEHDYF EFFSNDTSLY PHMYARAYIN FKNQEDIILF RDRFDGYVFL DNKGQEYPAI VEFAPFQKAA KKKTKKRDTK VGTIDDDPEY RKFLSYATD NEKMTSTPET LLEEIEAKNR ELIAKKTTP LLSFLKNKQRM REEKREERRR REIERKRQRE EERRKWKEEE KRKRKDIEKL KKIDRIPERD KLKDEPKIKV HR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com