

Product Data Sheet

Recombinant Human DOK7 Protein

Catalog #	Source	Reactivity	Applications
CST1815	E. coli	H	
Description	Recombinant Human DOK7 Protein is produced by E. coli expression system and the target gene encoding Val231-Ala380 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DOK7		
Alternative Names	C4orf25; Protein Dok-7; Downstream of tyrosine kinase 7		
Entrez Gene	285489		
SwissProt	Q18PE1		
Purity			
Chemical Structure	VAQEALETQ LEKRLSLLSH AGRPGSGGDD RSLSSSSSEA SHLDVSASSR LTAWPEQSSS SASTSQEGPR PAAAQAAGEA MVGASRPPPK PLRPRQLQEV GRQSSSDSGI ATGSHSSYSS SLSSYAGSSL DVWRATDELG SLLSLPAAGA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com