

Product Data Sheet

Recombinant Human SLITRK5 Protein

Catalog #	Source	Reactivity	Applications
CST1823	E. coli	H	
Description	Recombinant Human SLITRK5 Protein is produced by E. coli expression system and the target gene encoding Ser301-Tyr420 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SLITRK5		
Alternative Names	KIAA0918; LRRC11; SLIT and NTRK-like protein 5; Leucine-rich repeat-containing protein 11		
Entrez Gene	26050		
SwissProt	O94991		
Purity			
Chemical Structure	STTGYLHTTP ASVNSVATSS SAVYKPPLKP PKGTRQPNKP RVRPTSRQPS KDLGYSNYGP SIAYQTKSPV PLECPTACSC NLQISDLGLN VNCQERKIES IAELQPKPYN PKKMYLTENY		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com