

Product Data Sheet

Recombinant Human ENTPD3 Protein

Catalog #	Source	Reactivity	Applications
CST1824	E. coli	H	
Description	Recombinant Human ENTPD3 Protein is produced by E. coli expression system and the target gene encoding Leu51-Ser460 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ENTPD3		
Alternative Names	CD39L3; Ectonucleoside triphosphate diphosphohydrolase 3; NTPDase 3; CD39 antigen-like 3; Ecto-ATP diphosphohydrolase 3; Ecto-ATPDase 3; Ecto-ATPase 3; Ecto-apyrase 3; HB6		
Entrez Gene	956		
SwissProt	O75355		
Purity			
Chemical Structure	LPPGLKYGIV LDAGSSRTTV YVYQWPAEKE NNTGVVSQTF KCSVKGSGIS SYGNNPQDVP RAFEECMQKV KGQVPSHLHG STPIHLGATA GMRLRLQNE TAANEVLESI QSYFKSQPFD FRGAQIISGQ EEGVYGWITA NYLMGNFLEK NLWHMWVHPH GVETT GALDL GGASTQISFV AGEKMDLNTS DIMQVSPLYGY VYTLYTHSFQ CYGRNEAEKK FLAMLLQNSP TK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com