

Recombinant Human ADGRB1 Protein

Catalog #	Source	Reactivity	Applications
CST1826	E. coli	H	
Description	Recombinant Human ADGRB1 Protein is produced by E. coli expression system and the target gene encoding Gly1271-Thr1540 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	BAI1		
Alternative Names	Brain-specific angiogenesis inhibitor 1		
Entrez Gene	575		
SwissProt	O14514		
Purity			
Chemical Structure	GPPTNFNSLP ANVSKLHLHG SPRYPGGPLP DFPNHSLTK RDKAPKSSFV GDGDIFKKLD SELSRAQEKA LDTSYVILPT ATATLRPKPK EEPKYSIHID QMPQTRLIHL STAPEASLPA RSPPSRQPPS GGPPEAPPAQ PPPPPPPPPP PPQQPLPPPP NLEPAPPSLG DPGEPAAHPG PSTGPSTKNE NVATLSVSSL ERRKSRYAEL DFEKIMHTRK RHQDMFQDLN RK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com